

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti-NPS



Numéro de catalogue: 23609-1-AP

Informations de base

Numéro de catalogue:

23609-1-AP

Numéro d'acquisition GenBank:

BC156665

Méthode de purification:

Purifié par affinité contre l'antigène

Taille:

150ul, Concentration: 900 µg/ml by Nanodrop and 480 µg/ml by Bradford method using BSA as the standard;

Identification du gène (NCBI):

594857

Dilutions recommandées:

IHC 1:20-1:200

Hôte:

Lapin

Nom complet:

neuropeptide S

MW calculé

89 aa, 10 kDa

Isotype:

IgG

Immunogen Catalog Number:

AG20322

Applications

Applications testées:

IHC, ELISA

Spécificité de l'espèce:

Humain

Contrôles positifs:

IHC : tissu cérébral humain, tissu de cancer du pancréas humain, tissu pancréatique humain

Remarque-IHC: il est suggéré de démasquer l'antigène avec un tampon de **TE buffer pH 9.0; (*) À défaut, 'le démasquage de l'antigène peut être 'effectué avec un tampon citrate pH 6,0.**

Informations générales

Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20°C

*** Les 20ul contiennent 0,1% de BSA.

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com


W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Données de validation sélectionnées




Immunohistochemical analysis of paraffin-embedded human brain slide using 23609-1-AP (NPS Antibody) at dilution of 1:50 (under 10x lens).



Immunohistochemical analysis of paraffin-embedded human brain slide using 23609-1-AP (NPS Antibody) at dilution of 1:50 (under 40x lens).



Immunohistochemical analysis of paraffin-embedded human pancreas cancer tissue slide using 23609-1-AP (NPS antibody) at dilution of 1:100 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunohistochemical analysis of paraffin-embedded human pancreas cancer tissue slide using 23609-1-AP (NPS antibody) at dilution of 1:100 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).