

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti-C7orf63



Numéro de catalogue: 24001-1-AP

Informations de base

Numéro de catalogue: 24001-1-AP	Numéro d'acquisition GenBank: BC132771	Méthode de purification: Purifié par affinité contre l'antigène
Taille: 150ul , Concentration: 260 µg/ml by Nanodrop;	Identification du gène (NCBI): 79846	Dilutions recommandées: WB 1:500-1:1000
Hôte: Lapin	Nom complet: chromosome 7 open reading frame 63	
Isotype: IgG	MW calculé: 941 aa, 106 kDa	
Immunogen Catalog Number: AG21205	MW observés: 110 kDa	

Applications

Applications testées: WB, ELISA	Contrôles positifs: WB : cellules NIH/3T3,
Spécificité de l'espèce: souris	

Informations générales

Stockage

Stockage:
Stocker à -20°C. Stable pendant un an après l'expédition.
Tampon de stockage:
PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3
L'aliquotage n'est pas nécessaire pour le stockage à -20C

***** Les 20ul contiennent 0,1% de BSA.**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



NIH/3T3 cells were subjected to SDS PAGE followed by western blot with 24001-1-AP (C7orf63 antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.