

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti-FOXD4



Numéro de catalogue: 24835-1-AP

Informations de base

Numéro de catalogue:
24835-1-AP

Taille:
150ul, Concentration: 1000 µg/ml by
Nanodrop and 440 µg/ml by Bradford
method using BSA as the standard;

Hôte:
Lapin

Isotype:
IgG

Immunogen Catalog Number:
AG19226

Numéro d'acquisition GenBank:
BC136570

Identification du gène (NCBI):

Nom complet:
forkhead box D4

MW calculé
439 aa, 47 kDa

MW observés:
65-70 kDa

Méthode de purification:
Purification par affinité contre
l'antigène

Dilutions recommandées:
WB 1:500-1:2000

Applications

Applications testées:
WB, ELISA

Spécificité de l'espèce:
Humain, rat, souris

Contrôles positifs:

WB: cellules NIH/3T3, cellules Jurkat, cellules RAW
264.7, cellules U-937

Informations générales

Forkhead box protein D4 (FOXD4) also named as FKHL9 or Myeloid factor alpha is a 439 amino acid protein, which contains one fork head DNA-binding domain. FOXD4L6 localizes in the nucleus. FOXD4 as an embryonic transcriptional regulator involves in leukemogenesis. The molecular weight of FOXD4L6 is 47kDa, but our experiment always detected a 60-70kDa protein by western blot.

Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20°C

*** Les 20ul contiennent 0,1% de BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



NIH/3T3 cells were subjected to SDS PAGE followed by western blot with 24835-1-AP (FOXD4 Antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.