

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti-HOXA2



Numéro de catalogue: 25044-1-AP

Informations de base

Numéro de catalogue:

25044-1-AP

Taille:

150ul, Concentration: 600 µg/ml by Nanodrop and 360 µg/ml by Bradford method using BSA as the standard;

Hôte:

Lapin

Isotype:

IgG

Immunogen Catalog Number:

AG19224

Numéro d'acquisition GenBank:

BC130571

Identification du gène (NCBI):

3199

Nom complet:

homeobox A2

MW calculé

376 aa, 41 kDa

MW observés:

41-43 kDa

Méthode de purification:

Purification par affinité contre l'antigène

Dilutions recommandées:

WB 1:500-1:1000

Applications

Applications testées:

WB, ELISA

Spécificité de l'espèce:

Humain, rat

Contrôles positifs:

WB: cellules HepG2, cellules C6, cellules Jurkat

Informations générales

HOXA2, also named as Homeobox protein Hox-A2, is a 376 amino acid protein, which contains 1 homeobox DNA-binding domain and belongs to the Antp homeobox family. HOXA2 localizes in the nucleus and has an association with Microtia, hearing impairment, and cleft palate. HOXA2 as a sequence-specific transcription factor is part of a developmental regulatory system, which provides cells with specific positional identities on the anterior-posterior axis.

Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azote de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20°C

***** Les 20ul contiennent 0,1% de BSA.**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



HepG2 cells were subjected to SDS PAGE followed by western blot with 25044-1-AP (HOXA2 Antibody) at dilution of 1:600 incubated at 4 degree celsius over night.