

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti-FAM40B



Numéro de catalogue: 25163-1-AP

Phare

1 Publications

Informations de base

Numéro de catalogue:
25163-1-AP

Taille:
150ul, Concentration: 760 µg/ml by
Nanodrop and 260 µg/ml by Bradford
method using BSA as the standard;

Hôte:
Lapin

Isotype:
IgG

Immunogen Catalog Number:
AG18230

Numéro d'acquisition GenBank:
BC019064

Identification du gène (NCBI):
57464

Nom complet:
family with sequence similarity 40,
member B

MW calculé
834 aa, 95 kDa

MW observés:
95-100 kDa

Méthode de purification:
Purification par affinité contre
l'antigène

Dilutions recommandées:
WB 1:1000-1:4000

Applications

Applications testées:
WB, ELISA

Demandes citées:
IF, WB

Spécificité de l'espèce:
Humain, rat, souris

Espèces citées:
Humain

Contrôles positifs:

WB : tissu testiculaire de souris, tissu testiculaire de rat

Informations générales

FAM40B, also known as STRIP2, is a member of the striatin-interacting phosphatase and kinase (STRIPAK) complex that is involved in the regulation of various processes such as cell proliferation and differentiation. Fam40b in ESCs is as a perinuclear and nucleolar protein with a molecular weight of 96kDa (PMID: 25010986). The expression of Fam40b is essential for the lineage commitment of murine embryonic stem cells (mESCs) into differentiated somatic cells via mechanisms involving pluripotency and epigenetic networks. Catalog#17293-1-AP specially recognises the 96 kDa protein.

Publications notables

Autrice	Pubmed ID	Journal	Application
Li Dai	31093976	J Cell Physiol	WB,IF

Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20C

*** Les 20ul contiennent 0,1% de BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



Various lysates were subjected to SDS PAGE followed by western blot with 25163-1-AP (FAM40B antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.