

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti-NCOA6



Numéro de catalogue: 25241-1-AP

2 Publications

Informations de base

Numéro de catalogue:

25241-1-AP

Taille:

150ul, Concentration: 900 µg/ml by Nanodrop and 313 µg/ml by Bradford method using BSA as the standard;

Hôte:

Lapin

Isotype:

IgG

Immunogen Catalog Number:

AG18737

Numéro d'acquisition GenBank:

BC136272

Identification du gène (NCBI):

23054

Nom complet:

nuclear receptor coactivator 6

MW calculé

2063 aa, 219 kDa

MW observés:

68 kDa, 219 kDa, 110 kDa

Méthode de purification:

Purification par affinité contre l'antigène

Dilutions recommandées:

IP 0.5-4.0 ug for IP and 1:500-1:2000 for WB

Applications

Applications testées:

FC, IP, ELISA

Demandes citées:

IF, WB

Spécificité de l'espèce:

Humain, souris

Espèces citées:

Humain

Contrôles positifs:

IP: cellules HEK-293,

Informations générales

NCOA6 (Nuclear receptor coactivator 6), also named as AIB3 or ANTP, is a 2063 amino acid protein, which contains two LXXLL motifs. NCOA6 localizes in the nucleus and is highly expressed in brain, prostate, testis and ovary. NCOA6 directly binds nuclear receptors and stimulates the transcriptional activities in a hormone-dependent fashion. The molecular weight of NCOA6 is 219 kDa, but we always detected 68 kDa protein by western blot.

Publications notables

Autrice	Pubmed ID	Journal	Application
Yetao Xu	30195759	Mol Ther Nucleic Acids	WB
Liang Yan	28546420	Mol Pharmacol	WB, IF

Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20°C

*** Les 20ul contiennent 0,1% de BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



1X10⁶ HeLa cells were stained with 0.20ug NCOA6 antibody (25241-1-AP, red) and control antibody (blue). Fixed with 90% MeOH.



IP Result of anti-NCOA6 (IP:25241-1-AP, 4ug; Detection:25241-1-AP 1:1000) with HEK-293 cells lysate 2560ug.