

À des fins de recherche uniquement

# Anticorps Polyclonal de lapin anti-MRP1



Numéro de catalogue: 25340-1-AP

## Informations de base

Numéro de catalogue:  
25340-1-AP

Taille:  
150ul

Hôte:  
Lapin

Isotype:  
IgG

Numéro d'acquisition GenBank:  
NM\_004996

Identification du gène (NCBI):  
4363

Nom complet:  
ATP-binding cassette, sub-family C  
(CFTR/MRP), member 1

MW calculé  
173 kDa

Méthode de purification:  
Purification par affinité contre  
l'antigène

## Applications

Applications testées:  
ELISA

Spécificité de l'espèce:  
Humain

## Informations générales

Multidrug resistant-associated protein 1 (MRP1, also known as ABCC1) is a member of the ATP-binding cassette (ABC) transporter protein superfamily, subfamily C. MRP1 is overexpressed in cancers and contributes to the occurrence of multidrug resistance (MDR). Expression of MRP1 has been considered as a negative marker for the outcome of chemotherapy. Recently MRP1 has been reported to play a crucial role in A $\beta$  clearance at the blood-brain barrier.

## Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20C

\*\*\* Les 20ul contiennent 0,1% de BSA.

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

