

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti-SLC35F2



Numéro de catalogue: 25526-1-AP

Phare

3 Publications

Informations de base

Numéro de catalogue: 25526-1-AP	Numéro d'acquisition GenBank: BC039195	Méthode de purification: Purification par affinité contre l'antigène
Taille: 150ul , Concentration: 600 µg/ml by Nanodrop and 307 µg/ml by Bradford method using BSA as the standard;	Identification du gène (NCBI): 54733	Dilutions recommandées: WB 1:500-1:1000
Hôte: Lapin	Nom complet: solute carrier family 35, member F2	
Isotype: IgG	MW calculé: 374 aa, 41 kDa	
Immunogen Catalog Number: AG20901	MW observés: 41-50 kDa	

Applications

Applications testées: WB, ELISA	Contrôles positifs: WB : cellules A549, tissu cérébral de souris
Demandes citées: WB	
Spécificité de l'espèce: Humain, souris	
Espèces citées: Humain	

Informations générales

SLC35F2, also named as Solute carrier family 35 member F2, is a 374 amino acid protein, which belongs to the SLC35F solute transporter family. SLC35F2 as a Multi-pass membrane protein is a putative solute transporter. SLC35F2 is highly homologous to the lung squamous cell cancer-related gene, LSCC3, which is highly expressed in lung squamous cell tumour tissues. However, the clinical implication of the SLC35F2 gene in tumour development and progression remains unclear.

Publications notables

Autrice	Pubmed ID	Journal	Application
Arun Pandian Chandrasekaran	34815782	Theranostics	WB
Mei Chen	31849485	Onco Targets Ther	WB
Michael D Nyquist	28465296	Mol Cancer Res	WB

Stockage


Stockage:
Stocker à -20°C. Stable pendant un an après l'expédition.
Tampon de stockage:
PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3
L'aliquotage n'est pas nécessaire pour le stockage à -20C

*** Les 20ul contiennent 0,1% de BSA.

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



A549 cells were subjected to SDS PAGE followed by western blot with 25526-1-AP (SLC35F2 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.