

À des fins de recherche uniquement

# Anticorps Polyclonal de lapin anti-ECHDC2



Numéro de catalogue: 26126-1-AP

## Informations de base

Numéro de catalogue:

26126-1-AP

Taille:

150ul

Hôte:

Lapin

Isotype:

IgG

Immunogen Catalog Number:

AG23510

Numéro d'acquisition GenBank:

BC044574

Identification du gène (NCBI):

55268

Nom complet:

enoyl Coenzyme A hydratase domain containing 2

Méthode de purification:

Purification par affinité contre l'antigène

## Applications

Applications testées:

WB, ELISA

Spécificité de l'espèce:

Humain, souris

## Informations générales

### Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20C

**\*\*\* Les 20ul contiennent 0,1% de BSA.**

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Données de validation sélectionnées



mouse liver mitochondria tissue were subjected to SDS PAGE followed by western blot with 26126-1-AP (ECHDC2 Antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.