

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti-ZNF169



Numéro de catalogue: 26344-1-AP

Informations de base

Numéro de catalogue:

26344-1-AP

Taille:

150ul, Concentration: 450 µg/ml by Nanodrop and 313 µg/ml by Bradford method using BSA as the standard;

Hôte:

Lapin

Isotype:

IgG

Immunogen Catalog Number:

AG23805

Numéro d'acquisition GenBank:

BC047702

Identification du gène (NCBI):

169841

Nom complet:

zinc finger protein 169

MW observés:

68 kDa

Méthode de purification:

Purification par affinité contre l'antigène

Dilutions recommandées:

WB 1:500-1:1000

Applications

Applications testées:

WB, ELISA

Spécificité de l'espèce:

Humain, souris

Contrôles positifs:

WB : tissu d'intestin grêle de souris, cellules HeLa, cellules HepG2, tissu hépatique de souris

Informations générales

ZNF169 (Zinc finger protein 169) is a 603 amino acid nuclear protein that contains one KRAB domain and thirteen C2H2-type zinc fingers. ZNF169 is highly expressed in kidney and weakly expressed in spleen, liver, small intestine and heart, where it functions as a transcription regulator.

Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20C

*** Les 20ul contiennent 0,1% de BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



mouse small intestine tissue were subjected to SDS PAGE followed by western blot with 26344-1-AP (ZNF169 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.