

À des fins de recherche uniquement

# Anticorps Polyclonal de lapin anti-HOXA10



Numéro de catalogue: 26497-1-AP

Phare

5 Publications

## Informations de base

Numéro de catalogue:

26497-1-AP

Taille:

150ul, Concentration: 467 µg/ml by Nanodrop and 467 µg/ml by Bradford method using BSA as the standard;

Hôte:

Lapin

Isotype:

IgG

Immunogen Catalog Number:

AG24156

Numéro d'acquisition GenBank:

BC013971

Identification du gène (NCBI):

3206

Nom complet:

homeobox A10

MW calculé

393 aa, 41 kDa

MW observés:

45-50 kDa

Méthode de purification:

Purification par affinité contre l'antigène

Dilutions recommandées:

WB 1:500-1:1000

## Applications

Applications testées:

WB, ELISA

Demandes citées:

IF, IHC, WB

Spécificité de l'espèce:

Humain

Espèces citées:

bovin, Humain, singe

Contrôles positifs:

WB : cellules HEK-293, cellules HeLa

## Informations générales

HOXA10, also named as Homeobox protein Hox-A10, is a 410 amino acid protein, which contains 1 homeobox DNA-binding domain and belongs to the Abd-B homeobox family. HOXA10 localizes in the nucleus and can interact with SIRT2. HOXA10 as sequence-specific transcription factor, is part of a developmental regulatory system that provides cells with specific positional identities on the anterior-posterior axis. HoxA10 expression increases during the midsecretory phase of the menstrual cycle, which corresponds with increased levels of circulating progesterone.

## Publications notables

Autrice	Pubmed ID	Journal	Application
Jing Guo	32183826	Reprod Biol Endocrinol	WB, IHC
Peipei Yan	28116645	J Assist Reprod Genet	WB
Shuzhen Guan	36562948	Front Med	WB

## Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azotate de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20°C

\*\*\* Les 20ul contiennent 0,1% de BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Données de validation sélectionnées



HEK-293 cells were subjected to SDS PAGE followed by western blot with 26497-1-AP (HOXA10 Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.