

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti-MUC7



Numéro de catalogue: 26544-1-AP

Informations de base

Numéro de catalogue:
26544-1-AP

Taille:
150ul, Concentration: 300 µg/ml by
Nanodrop;

Hôte:
Lapin

Isotype:
IgG

Immunogen Catalog Number:
AG24216

Numéro d'acquisition GenBank:
BC025688

Identification du gène (NCBI):
4589

Nom complet:
mucin 7, secreted

MW calculé
39 kDa

MW observés:
38-40 kDa

Méthode de purification:
Purification par affinité contre
l'antigène

Dilutions recommandées:
WB 1:500-1:2000

Applications

Applications testées:
WB, ELISA
Spécificité de l'espèce:
Humain

Contrôles positifs:
WB : salive humain,

Informations générales

Stockage


Stockage:
Stocker à -20°C. Stable pendant un an après l'expédition.
Tampon de stockage:
PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3
L'aliquotage n'est pas nécessaire pour le stockage à -20C

***** Les 20ul contiennent 0,1% de BSA.**

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA) E: proteintech@ptglab.com W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



human saliva were subjected to SDS PAGE followed by western blot with 26544-1-AP (MUC7 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.