

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti-GNG3



Numéro de catalogue: 26892-1-AP

Informations de base

Numéro de catalogue:	26892-1-AP	Numéro d'acquisition GenBank:	BC015563	Méthode de purification:	Purification par affinité contre l'antigène
Taille:	150ul, Concentration: 600 µg/ml by Nanodrop;	Identification du gène (NCBI):	2785	Dilutions recommandées:	WB 1:500-1:3000
Hôte:	Lapin	Nom complet:	guanine nucleotide binding protein (G protein), gamma 3		
Isotype:	IgG	MW calculé	8 kDa		
Immunogen Catalog Number:	AG25341	MW observés:	10 kDa		

Applications

Applications testées:	WB, ELISA	Contrôles positifs:	WB : tissu cérébral de souris, tissu cérébral de rat
Spécificité de l'espèce:	rat, souris		

Informations générales

GNG3, also named as GNGT3, belongs to the G protein gamma family. G proteins are heterotrimers of alpha, beta, and gamma subunits. Gamma subunits, such as GNG3, contribute to the specificity of the hundreds of receptor signaling pathways involving G proteins. GNG3 is predominantly expressed in the brain, where its expression increases during postnatal development.

Stockage


Stockage:
Stocker à -20°C. Stable pendant un an après l'expédition.
Tampon de stockage:
PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3
L'aliquotage n'est pas nécessaire pour le stockage à -20C

*** Les 20ul contiennent 0,1% de BSA.

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



Various lysates were subjected to SDS PAGE followed by western blot with 26892-1-AP (GNG3 antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.