

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti-FC tag



Numéro de catalogue: 28108-1-AP

Informations de base

Numéro de catalogue:

28108-1-AP

Taille:

150ul, Concentration: 150 µg/ml by Nanodrop;

Hôte:

Lapin

Isotype:

IgG

Immunogen Catalog Number:

AG27844

Numéro d'acquisition GenBank:

Identification du gène (NCBI):

Nom complet:

MW calculé

26 kDa

MW observés:

55 kDa

Méthode de purification:

Purification par affinité contre l'antigène

Dilutions recommandées:

WB 1:2000-1:10000

Applications

Applications testées:

WB, ELISA

Spécificité de l'espèce:

souris

Contrôles positifs:

WB : Protéine recombinante,

Informations générales

This antibody detects the Fc fragment of mouse IgG2a.

Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20C

*** Les 20ul contiennent 0,1% de BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



Various lysates were subjected to SDS PAGE followed by western blot with 28108-1-AP (FC tag antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.