

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti-CDKN2A/p19ARF



Numéro de catalogue: 28416-1-AP

Informations de base

Numéro de catalogue:

28416-1-AP

Taille:

150ul, Concentration: 500 µg/ml by Nanodrop;

Hôte:

Lapin

Isotype:

IgG

Immunogen Catalog Number:

AG28927

Numéro d'acquisition GenBank:

NM_001040654

Identification du gène (NCBI):

12578

Nom complet:

cyclin-dependent kinase inhibitor 2A

MW calculé

18 kDa

MW observés:

19 kDa

Méthode de purification:

Purification par affinité contre l'antigène

Dilutions recommandées:

WB 1:500-1:2000

Applications

Applications testées:

WB, ELISA

Spécificité de l'espèce:

souris

Contrôles positifs:

WB : cellules NIH/3T3, cellules C6, cellules Neuro-2a, cellules RAW 264.7

Informations générales

Cdkn2a (cyclin-dependent kinase inhibitor 2a) is a cell cycle regulating gene encoding alternate proteins, p16Ink4a and p19Arf, that arrest the G1 phase through CDK4/6 inhibition. As they are tumor suppressors, p16Ink4a and p19Arf can be targeted for treatment of various cancers, including melanoma, cervical, and esophageal cancers.

Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20C

*** Les 20ul contiennent 0,1% de BSA.

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



Various lysates were subjected to SDS PAGE followed by western blot with 28416-1-AP (Cdkn2a antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.