

À des fins de recherche uniquement

# Anticorps Polyclonal de lapin anti-TRPS1



Numéro de catalogue: 28753-1-AP **Phare**

## Informations de base

|   |   |  |
|---|---|--|
| <b>Numéro de catalogue:</b><br>28753-1-AP                       | <b>Numéro d'acquisition GenBank:</b><br>NM_014112       | <b>Méthode de purification:</b><br>Purification par affinité contre l'antigène |
| <b>Taille:</b><br>150ul , Concentration: 600 µg/ml by Nanodrop; | <b>Identification du gène (NCBI):</b><br>7227           | <b>Dilutions recommandées:</b><br>WB 1:500-1:3000                              |
| <b>Hôte:</b><br>Lapin   | <b>Nom complet:</b><br>trichorhinophalangeal syndrome I |  |
| <b>Isotype:</b><br>IgG  | <b>MW calculé:</b><br>142 kDa                           |  |
| <b>Immunogen Catalog Number:</b><br>AG30188                     | <b>MW observés:</b><br>142 kDa                          |  |

## Applications

|   |  |
|---|--|
| <b>Applications testées:</b><br>WB, ELISA | <b>Contrôles positifs:</b><br>WB : cellules HeLa, cellules MCF-7 |
| <b>Spécificité de l'espèce:</b><br>Humain |  |

## Informations générales

TRPS1 is a zinc finger transcriptional repressor involved in the regulation of chondrocyte and perichondrium development, containing a GATA-type zinc finger through which it binds to DNA. with nine zinc-finger domains and two C-terminal Ikaros-like zinc fingers. Its repressible function was dependent on the integrity of the Trps1 GATA-type zinc-finger domain and also required the C-terminal 119 amino acids of the protein, which harbor the two Ikaros-like zinc-finger domains

## Stockage

**Stockage:**  
Stocker à -20°C. Stable pendant un an après l'expédition.  
**Tampon de stockage:**  
PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3  
L'aliquotage n'est pas nécessaire pour le stockage à -20C

**\*\*\* Les 20ul contiennent 0,1% de BSA.**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Données de validation sélectionnées



Various lysates were subjected to SDS PAGE followed by western blot with 28753-1-AP (TRPS1 antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.



WB result of TRPS1 antibody (28753-1-AP; 1:2000; incubated at room temperature for 1.5 hours) with sh-Control and sh-TRPS1 transfected HeLa cells.