

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti-Phospho-MAP3K1 (Thr1400)



Numéro de catalogue: 28844-1-AP

1 Publications

Informations de base

Numéro de catalogue:

28844-1-AP

Taille:

100ul, Concentration: 300 µg/ml by Nanodrop;

Hôte:

Lapin

Isotype:

IgG

Numéro d'acquisition GenBank:

NM_005921

Identification du gène (NCBI):

4214

Nom complet:

mitogen-activated protein kinase kinase kinase 1

MW calculé

164 kDa

MW observés:

164 kDa

Méthode de purification:

Purification par affinité contre l'antigène

Dilutions recommandées:

WB 1:2000-1:10000

Applications

Applications testées:

WB, ELISA

Demandes citées:

WB

Spécificité de l'espèce:

Humain

Espèces citées:

rat

Contrôles positifs:

WB : cellules HeLa traitées à la calyculine A,

Informations générales

MAP3K1, also named as MAPKKK1, MEKK and MEKK1, belongs to the protein kinase superfamily, STE Ser/Thr protein kinase family and MAP kinase kinase kinase subfamily. MAP3K1 is a component of a protein kinase signal transduction cascade. It activates CHUK and IKBKB, the central protein kinases of the NF-kappa-B pathway. MAP3K1 activates the ERK and JNK kinase pathways by phosphorylation of MAP2K1 and MAP2K4. MAP3K1 is activated by autophosphorylation on Thr-1400 and Thr-1412(Uniprot).

Publications notables

Autrice	Pubmed ID	Journal	Application
Song Liu	37149277	Eur J Pharmacol	WB

Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20C

*** Les 20ul contiennent 0,1% de BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



Non-treated and Calyculin A treated HeLa cells were subjected to SDS PAGE followed by western blot with 28844-1-AP (Phospho-MAP3K1 (Thr1400) antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.