

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti-RIPK3



Numéro de catalogue: 29080-1-AP

Informations de base

Numéro de catalogue: 29080-1-AP	Numéro d'acquisition GenBank: BC062584	Méthode de purification: Purification par affinité contre l'antigène
Taille: 150ul , Concentration: 800 µg/ml by Nanodrop;	Identification du gène (NCBI): 11035	Dilutions recommandées: WB 1:2000-1:10000
Hôte: Lapin	Nom complet: receptor-interacting serine-threonine kinase 3	
Isotype: IgG	MW calculé 518 aa, 57 kDa	
Immunogen Catalog Number: AG30337	MW observés: 57 kDa	

Applications

Applications testées: WB, ELISA	Contrôles positifs: WB : cellules THP-1,
Spécificité de l'espèce: Humain	

Informations générales

RIPK3, also named as RIP3, a Ser/Thr kinase of RIP (Receptor Interacting Protein) family, is recruited to the TNFR1 signaling complex through RIP and has been shown to mediate apoptosis induction and NF-κB activation. RIPK3 is a nucleocytoplasmic shuttling protein and its unconventional nuclear localization signal (NLS, 442-472 aa) is sufficient to trigger apoptosis in the nucleus (PMID:18533105). It has 3 isoforms produced by alternative splicing. RIPK3 might form a homodimer within cells, and its apoptotic activity may not be required for this dimerization (PMID:18533105).

Stockage

Stockage:
Stocker à -20°C. Stable pendant un an après l'expédition.
Tampon de stockage:
PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3
L'aliquotage n'est pas nécessaire pour le stockage à -20C

*** Les 20ul contiennent 0,1% de BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



THP-1 cells were subjected to SDS PAGE followed by western blot with 29080-1-AP (RIPK3 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.