

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti- Phospho-S6 Ribosomal protein (Ser240/244)



Numéro de catalogue: **29220-1-AP**

Informations de base

Numéro de catalogue: 29220-1-AP	Numéro d'acquisition GenBank: NM_001010	Méthode de purification: Purification par affinité contre l'antigène
Taille: 100ul , Concentration: 800 µg/ml by Nanodrop;	Identification du gène (NCBI): 6194	Dilutions recommandées: WB 1:1000-1:6000
Hôte: Lapin	Nom complet: ribosomal protein S6	
Isotype: IgG	MW observés: 32 kDa	

Applications

Applications testées:
WB, ELISA
Spécificité de l'espèce:
Humain, souris

Contrôles positifs:
WB : cellules NIH/3T3 traitées au sérum foetal de bovin
à 10 %

Informations générales

Ribosomal protein S6 (RPS6), also named as Phosphoprotein NP33, may play an important role in controlling cell growth and proliferation through the selective translation of particular classes of mRNA. Ribosomal protein S6 is the major substrate of protein kinases in eukaryote ribosomes. The phosphorylation is stimulated by growth factors, tumor promoting agents, and mitogens and is dephosphorylated at growth arrest.

Stockage


Stockage:
Stocker à -20°C. Stable pendant un an après l'expédition.
Tampon de stockage:
PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3
L'aliquotage n'est pas nécessaire pour le stockage à -20C

***** Les 20ul contiennent 0,1% de BSA.**

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



Serum-free treated NIH/3T3 cells and 10 % fetal bovine serum treated NIH/3T3 cells were subjected to SDS PAGE followed by western blot with 29220-1-AP (Phospho-S6 Ribosomal protein (Ser240/244) antibody) at dilution of 1:3000 incubated at 4°C overnight. The membrane was stripped and re-blotted with GAPDH antibody as loading control.