

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti-WDR47



Numéro de catalogue: 29364-1-AP

Informations de base

Numéro de catalogue:

29364-1-AP

Taille:

150ul , Concentration: 900 µg/ml by Nanodrop;

Hôte:

Lapin

Isotype:

IgG

Immunogen Catalog Number:

AG30853

Numéro d'acquisition GenBank:

BC039254

Identification du gène (NCBI):

22911

Nom complet:

WD repeat domain 47

MW calculé

919 aa, 102 kDa

MW observés:

100 kDa

Méthode de purification:

Purification par affinité contre l'antigène

Dilutions recommandées:

WB 1:2000-1:16000

Applications

Applications testées:

WB, ELISA

Spécificité de l'espèce:

Humain

Contrôles positifs:

WB : cellules HeLa, cellules K-562

Informations générales

Stockage

Stockage:

Stocker à -20 °C

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20C

***** Les 20ul contiennent 0,1% de BSA.**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



Various lysates were subjected to SDS PAGE followed by western blot with 29364-1-AP (WDR47 antibody) at dilution of 1:8000 incubated at room temperature for 1.5 hours.