

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti-GEN1



Numéro de catalogue: 29617-1-AP

Informations de base

Numéro de catalogue:

29617-1-AP

Numéro d'acquisition GenBank:

NM_182625

Méthode de purification:

Purification par affinité contre l'antigène

Taille:

150ul, Concentration: 400 µg/ml by Nanodrop;

Identification du gène (NCBI):

348654

Dilutions recommandées:

WB 1:1000-1:4000

Hôte:

Lapin

Nom complet:

Gen homolog 1, endonuclease (Drosophila)

Isotype:

IgG

MW calculé

103 kDa

Immunogen Catalog Number:

AG31227

MW observés:

110 kDa

Applications

Applications testées:

WB, ELISA

Contrôles positifs:

WB : cellules HL-60, cellules MCF-7

Spécificité de l'espèce:

Humain

Informations générales

GEN1 is a endonuclease which resolves Holliday junctions (HJs) by the introduction of symmetrically related cuts across the junction point, to produce nicked duplex products in which the nicks can be readily ligated. GEN1 efficiently cleaves both single and double HJs contained within large recombination intermediates.

Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20C

*** Les 20ul contiennent 0,1% de BSA.

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



Various lysates were subjected to SDS PAGE followed by western blot with 29617-1-AP (GEN1 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.