

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti-SRGAP3



Numéro de catalogue: 29624-1-AP

Informations de base

Numéro de catalogue: 29624-1-AP	Numéro d'acquisition GenBank: NM_014850	Méthode de purification: Purification par affinité contre l'antigène
Taille: 150ul , Concentration: 550 µg/ml by Nanodrop;	Identification du gène (NCBI): 9901	Dilutions recommandées: WB 1:2000-1:12000 IF 1:50-1:500
Hôte: Lapin	Nom complet: SLIT-ROBO Rho GTPase activating protein 3	
Isotype: IgG	MW calculé: 124 kDa	
Immunogen Catalog Number: AG31196	MW observés: 125 kDa	

Applications

Applications testées: IF, WB, ELISA	Contrôles positifs:
Spécificité de l'espèce: Humain, rat, souris	WB : tissu cérébral de souris, tissu cérébral de rat IF : tissu cérébral de souris,

Informations générales

SRGAP3 (SLIT-ROBO Rho GTPase Activating Protein 3), also named MEGAP and WRP, is a member of the Slit-Robo sub-family of Rho GTPase-activating proteins (Rho GAPs), controls actin and microtubule dynamics through negative regulation of Rac (PMID: 23108406). Diseases associated with SRGAP3 include Pilomyxoid Astrocytoma and Avoidant Personality Disorder. SRGAP3 is highly expressed in fetal and adult brain, including the cortex and the hippocampus.

Stockage

Stockage:
Stocker à -20 °C.
Tampon de stockage:
PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3
L'aliquotage n'est pas nécessaire pour le stockage à -20C

***** Les 20ul contiennent 0,1% de BSA.**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



Various lysates were subjected to SDS PAGE followed by western blot with 29624-1-AP (SRGAP3 antibody) at dilution of 1:6000 incubated at room temperature for 1.5 hours.



Immunofluorescent analysis of (4% PFA) fixed mouse brain tissue using SRGAP3 antibody (29624-1-AP) at dilution of 1:200 and CoraLite® 488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).