

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti-KCNK12



Numéro de catalogue: 29762-1-AP

Informations de base

Numéro de catalogue:

29762-1-AP

Taille:

150ul, Concentration: 500 µg/ml by Nanodrop;

Hôte:

Lapin

Isotype:

IgG

Immunogen Catalog Number:

AG31561

Numéro d'acquisition GenBank:

NM_022055

Identification du gène (NCBI):

56660

Nom complet:

potassium channel, subfamily K, member 12

MW calculé

47KD

MW observés:

49 kDa

Méthode de purification:

Purification par affinité contre l'antigène

Dilutions recommandées:

WB 1:500-1:2000

Applications

Applications testées:

WB, ELISA

Spécificité de l'espèce:

Humain, souris

Contrôles positifs:

WB : tissu rénal de souris,

Informations générales

Potassium channel, subfamily K, member 12, also known as KCNK12, is a potassium channel containing two pore-forming P domains. The product of this gene has not been shown to be a functional channel, however, it may require other non-pore-forming proteins for activity. The calculated MW is 49 kDa, 29762-1-AP can detect a band around 49 kDa.

Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20C

*** Les 20ul contiennent 0,1% de BSA.

For technical support and original validation data for this product please contact:


T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com

W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



mouse kidney tissue were subjected to SDS PAGE followed by western blot with 29762-1-AP (KCNK12 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.