

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti- CDT1



Numéro de catalogue: 29766-1-AP

Informations de base

Numéro de catalogue:
29766-1-AP

Taille:
150ul , Concentration: 450 µg/ml by
Nanodrop;

Hôte:
Lapin

Isotype:
IgG

Immunogen Catalog Number:
AG31758

Numéro d'acquisition GenBank:
NM_030928

Identification du gène (NCBI):
81620

Nom complet:
chromatin licensing and DNA
replication factor 1

MW calculé
60KD

MW observés:
65-70 kDa

Méthode de purification:
Purification par affinité contre
l'antigène

Dilutions recommandées:
WB 1:1000-1:4000

Applications

Applications testées:
WB, ELISA

Spécificité de l'espèce:
Humain

Contrôles positifs:

WB : cellules MCF-7, cellules HeLa, cellules HepG2

Informations générales

Stockage

Stockage:

Stocker à -20 °C

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20C

***** Les 20ul contiennent 0,1% de BSA.**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



Various lysates were subjected to SDS PAGE followed by western blot with 29766-1-AP (CDT1 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.