

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti-ZNF643



Numéro de catalogue: 30018-1-AP

Informations de base

Numéro de catalogue:

30018-1-AP

Taille:

150ul, Concentration: 550 µg/ml by Nanodrop;

Hôte:

Lapin

Isotype:

IgG

Immunogen Catalog Number:

AG32346

Numéro d'acquisition GenBank:

BC017498

Identification du gène (NCBI):

65243

Nom complet:

zinc finger protein 643

MW observés:

70-75 kDa

Méthode de purification:

Purification par affinité contre l'antigène

Dilutions recommandées:

WB 1:500-1:3000

Applications

Applications testées:

WB, ELISA

Spécificité de l'espèce:

Humain

Contrôles positifs:

WB : cellules HepG2, cellules HuH-7, cellules SMMC-7721, cellules THP-1

Informations générales

ZNF643 (Zinc finger protein 643) is also named as ZFP69B. ZNF643, a putative transcription factor gene, is situated next to ZFP69, which has been linked to pathogenesis of human diabetes, as its allelic variation associates with impaired lipid storage in white adipose tissue (PMID: 19578398). It may be involved in transcriptional regulation, and is essential for Golgi structural integrity (PMID:29851555). ZNF643 has two isoforms with molecular weights of 50 kDa and 62 kDa. ZNF643 is modified by SUMO.

Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azote de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20°C

*** Les 20ul contiennent 0,1% de BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



Various lysates were subjected to SDS PAGE followed by western blot with 30018-1-AP (ZNF643 antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.