

À des fins de recherche uniquement

# Anticorps Polyclonal de lapin anti-GSDMA



Numéro de catalogue: 30354-1-AP

## Informations de base

|                                              |                                |                                             |
|----------------------------------------------|--------------------------------|---------------------------------------------|
| Numéro de catalogue:                         | Numéro d'acquisition GenBank:  | Méthode de purification:                    |
| 30354-1-AP                                   | BC109197                       | Purification par affinité contre l'antigène |
| Taille:                                      | Identification du gène (NCBI): | Dilutions recommandées:                     |
| 150ul, Concentration: 550 µg/ml by Nanodrop; | 284110                         | WB 1:2000-1:14000                           |
| Hôte:                                        | Nom complet:                   |                                             |
| Lapin                                        | gasdermin A                    |                                             |
| Isotype:                                     | MW calculé                     |                                             |
| IgG                                          | 445 aa, 49 kDa                 |                                             |
| Immunogen Catalog Number:                    | MW observés:                   |                                             |
| AG31284                                      | 49 kDa                         |                                             |

## Applications

|                          |                                                                                           |
|--------------------------|-------------------------------------------------------------------------------------------|
| Applications testées:    | Contrôles positifs:                                                                       |
| WB, ELISA                | WB : cellules COLO 320, cellules Daudi, tissu cutané de souris, tissu d'estomac de souris |
| Spécificité de l'espèce: |                                                                                           |
| Humain, souris           |                                                                                           |

## Informations générales

Gasdermins, a family of five pore-forming proteins (GSDMA-GSDME) in humans expressed predominantly in the skin, mucosa and immune sentinel cells, are key executioners of inflammatory cell death (pyroptosis), which recruits immune cells to infection sites and promotes protective immunity. GSDMA was predicted to enable phosphatidylinositol-4,5-bisphosphate binding activity, phosphatidylinositol-4-phosphate binding activity, and phosphatidylserine binding activity. Involved in apoptotic process.

## Stockage


Stockage:  
Stocker à -20°C. Stable pendant un an après l'expédition.  
Tampon de stockage:  
PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3  
L'aliquotage n'est pas nécessaire pour le stockage à -20C

\*\*\* Les 20ul contiennent 0,1% de BSA.

For technical support and original validation data for this product please contact:  
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)  
E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Données de validation sélectionnées



Various lysates were subjected to SDS PAGE followed by western blot with 30354-1-AP (GSDMA antibody) at dilution of 1:7000 incubated at room temperature for 1.5 hours.