

À des fins de recherche uniquement

# Anticorps Polyclonal de lapin anti-TREM1



Numéro de catalogue: 30496-1-AP

## Informations de base

<b>Numéro de catalogue:</b> 30496-1-AP	<b>Numéro d'acquisition GenBank:</b> BC017773	<b>Méthode de purification:</b> Purification par affinité contre l'antigène
<b>Taille:</b> 150ul , Concentration: 350 µg/ml by Nanodrop;	<b>Identification du gène (NCBI):</b> 54210	<b>Dilutions recommandées:</b> WB 1:500-1:1000
<b>Hôte:</b> Lapin	<b>Nom complet:</b> triggering receptor expressed on myeloid cells 1	
<b>Isotype:</b> IgG	<b>MW calculé:</b> 234 aa, 26 kDa	
<b>Immunogen Catalog Number:</b> AG32967	<b>MW observés:</b> 30 kDa	

## Applications

<b>Applications testées:</b> WB, ELISA	<b>Contrôles positifs:</b> WB : cellules THP-1,
<b>Spécificité de l'espèce:</b> Humain	

## Informations générales

Triggering receptor expressed on myeloid cells 1 (TREM1, also known as CD354) is a member of the immunoglobulin superfamily, which is expressed on neutrophils and monocytes. TREM1 plays a role in the innate immune system. TREM1 amplifies neutrophil and monocyte-mediated inflammatory responses triggered by bacterial and fungal infections by stimulating release of pro-inflammatory chemokines and cytokines, as well as increased surface expression of cell activation markers.

## Stockage

**Stockage:**  
Stocker à -20°C. Stable pendant un an après l'expédition.  
**Tampon de stockage:**  
PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3  
L'aliquotage n'est pas nécessaire pour le stockage à -20C

**\*\*\* Les 20ul contiennent 0,1% de BSA.**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Données de validation sélectionnées



THP-1 cells were subjected to SDS PAGE followed by western blot with 30496-1-AP (TREM1 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.