

À des fins de recherche uniquement

# Anticorps Polyclonal de lapin anti-CRISP3



Numéro de catalogue: 51022-1-AP

## Informations de base

Numéro de catalogue:

51022-1-AP

Taille:

150ul, Concentration: 133 µg/ml by Bradford method using BSA as the standard;

Hôte:

Lapin

Isotype:

IgG

Numéro d'acquisition GenBank:

BC069602

Identification du gène (NCBI):

10321

Nom complet:

cysteine-rich secretory protein 3

MW calculé

28 kDa

Méthode de purification:

Purification par affinité contre l'antigène

## Applications

Applications testées:

WB, ELISA

Spécificité de l'espèce:

Humain

## Informations générales

CRISP3, also known as SGP28, is a glycoprotein that belongs to a family of cysteine-rich secretory proteins (CRISPs). CRISP3 was originally discovered in human neutrophils, but transcripts are widely distributed in exocrine glands (salivary glands, pancreas, and prostate) and also found at lower levels in epididymis, ovary, thymus, and colon. CRISP3 is a potential biomarker for prostate cancer.

## Stockage

Stockage:

Stocker à -20 °C

Tampon de stockage:

PBS avec azote de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20C

\*\*\* Les 20ul contiennent 0,1% de BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Données de validation sélectionnées



HEK-293 cells were subjected to SDS PAGE followed by western blot with 51022-1-AP (CRISP3 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.