

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti-BNC1



Numéro de catalogue: 55233-1-AP

Informations de base

Numéro de catalogue:

55233-1-AP

Taille:

150ul , Concentration: 413 µg/ml by Bradford method using BSA as the standard;

Hôte:

Lapin

Isotype:

IgG

Numéro d'acquisition GenBank:

NM_001717

Identification du gène (NCBI):

646

Nom complet:

basonuclin 1

MW calculé

111 kDa

Méthode de purification:

Purification par affinité contre l'antigène

Dilutions recommandées:

IF 1:50-1:500

Applications

Applications testées:

IF, ELISA

Spécificité de l'espèce:

Humain, souris

Contrôles positifs:

IF : cellules HepG2,

Informations générales

Stockage

Stockage:

Stocker à -20 °C

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20C

*** Les 20ul contiennent 0,1% de BSA.

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



Immunofluorescent analysis of (4% PFA) fixed HepG2 cells using 55233-1-AP (BNC1 antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).