

À des fins de recherche uniquement

Anticorps Polyclonal de lapin anti-MUC13



Numéro de catalogue: 55487-1-AP

Informations de base

Numéro de catalogue:

55487-1-AP

Taille:

150ul, Concentration: 1500 µg/ml by 56667

Nanodrop;

Hôte:

Lapin

Isotype:

IgG

Numéro d'acquisition GenBank:

NM_033049

Identification du gène (NCBI):

56667

Nom complet:

mucin 13, cell surface associated

MW calculé

55 kDa

MW observés:

90 kDa

Méthode de purification:

Purification par affinité contre l'antigène

Dilutions recommandées:

WB 1:500-1:2000

Applications

Applications testées:

WB, ELISA

Spécificité de l'espèce:

Humain

Contrôles positifs:

WB : cellules HCT 116, cellules HEK-293

Informations générales

Mucins (MUC) are highly glycosylated proteins that are abundantly expressed in mammalian epithelial cells. Mucin 13 (MUC13) is a member of the transmembrane glycoprotein mucins which is normally expressed in the large intestine, trachea, kidney, small intestine, and gastric epithelium (PMID: 22027689). MUC13 is involved in the tumorigenesis of multiple malignancies, including pancreatic cancer, colorectal cancer, renal cancer and lung cancer (PMID: 22027689; PMID: 31427737; PMID: 28205224; PMID: 34966453).

Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20°C

*** Les 20ul contiennent 0,1% de BSA


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



Various lysates were subjected to SDS PAGE followed by western blot with 55487-1-AP (MUC13 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.