

À des fins de recherche uniquement

Anticorps Monoclonal anti-NUMBL

Numéro de catalogue: 60009-1-Ig **2 Publications**



Informations de base

Numéro de catalogue: 60009-1-Ig	Numéro d'acquisition GenBank: BC001794	Méthode de purification: Précipitation de l'acide caprylique/du sulfate d'ammonium
Taille: 150ul	Identification du gène (NCBI): 9253	CloneNo.: 4C10C5
Hôte: Mouse	Nom complet: numb homolog (Drosophila)-like	Dilutions recommandées: WB 1:200-1:1000
Isotype: IgG2a	MW calculé 65 kDa	
Immunogen Catalog Number: AG0156		

Applications

Applications testées: WB, ELISA	Contrôles positifs: WB : cellules HepG2,
Demandes citées: WB	
Spécificité de l'espèce: Humain	
Espèces citées: Humain	

Informations générales

Numb-like (Numbl or Nbl) binds Notch and is an intracellular modulator of Notch signaling pathway. Numbl contains a phosphotyrosine binding (PTB) domain that could bind the Y682 motif present in the cytoplasmic tail of APP. Deduced molecular weight of Numbl is 72 kDa.

Publications notables

Autrice	Pubmed ID	Journal	Application
Xingyu Du	36335371	J Transl Med	WB
Shu-Ting Pan	34484337	J Oncol	WB

Stockage


Stockage:
Stocker à -20°C. Stable pendant un an après l'expédition.
Tampon de stockage:
PBS avec azoture de sodium à 0,1 % et glycérol à 50 % pH 7,3
L'aliquotage n'est pas nécessaire pour le stockage à -20C

*** Les 20ul contiennent 0,1% de BSA.

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



HepG2 cells were subjected to SDS PAGE followed by western blot with 60009-1-Ig (NUMBL antibody) at dilution of 1:300 incubated at room temperature for 1.5 hours.