

À des fins de recherche uniquement

# Anticorps Monoclonal anti-Inhibin beta A



Numéro de catalogue:60015-2-Ig

## Informations de base

Numéro de catalogue: 60015-2-Ig	Numéro d'acquisition GenBank: BC007858	Méthode de purification: Précipitation de l'acide caprylique/du sulfate d'ammonium
Taille: 150ul , Concentration: 1000 µg/ml by Bradford method using BSA as the standard;	Identification du gène (NCBI): 3624	CloneNo.: 1E10C7
Hôte: Mouse	Nom complet: inhibin, beta A	
Isotype: IgG1	MW calculé 47 kDa	
Immunogen Catalog Number: AG1001		

## Applications

Applications testées:  
ELISA  
Spécificité de l'espèce:  
Humain

## Informations générales

Inhibin  $\beta$ A (INHBA) is a ligand in the transforming growth factor  $\beta$  (TGF- $\beta$ ) superfamily. Inhibin  $\beta$ A forms a disulfide-linked homo-dimer known as activin A, which has been implicated in multiple biologic processes, including neoplastic progression. The inhibin  $\beta$ A can also bind the  $\beta$ B and  $\alpha$  isoforms to form activin AB and inhibin, respectively.

## Stockage

Stockage:  
Stocker à -20°C. Stable pendant un an après l'expédition.  
Tampon de stockage:  
PBS avec azoture de sodium à 0,1 % et glycérol à 50 % pH 7,3  
L'aliquotage n'est pas nécessaire pour le stockage à -20C

**\*\*\* Les 20ul contiennent 0,1% de BSA.**

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

