

À des fins de recherche uniquement

# Anticorps Monoclonal anti-TWEAKR

Numéro de catalogue: 60104-1-Ig



## Informations de base

Numéro de catalogue:

60104-1-Ig

Taille:

150ul

Hôte:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG2689

Numéro d'acquisition GenBank:

BC002718

Identification du gène (NCBI):

51330

Nom complet:

tumor necrosis factor receptor  
superfamily, member 12A

MW calculé

129 aa, 14 kDa

Méthode de purification:

Précipitation de l'acide caprylique/du  
sulfate d'ammonium

CloneNo.:

2E10C10

## Applications

Applications testées:

WB, ELISA

Spécificité de l'espèce:

Humain

## Informations générales

### Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azoture de sodium à 0,1 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20C

**\*\*\* Les 20ul contiennent 0,1% de BSA.**


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Données de validation sélectionnées



human heart tissue were subjected to SDS PAGE followed by western blot with 60104-1-Ig (TWEAKR antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.