

À des fins de recherche uniquement

# Anticorps Monoclonal anti-BTBD10

Numéro de catalogue: 60129-1-Ig



## Informations de base

<b>Numéro de catalogue:</b> 60129-1-Ig	<b>Numéro d'acquisition GenBank:</b> BC017072	<b>Méthode de purification:</b> Précipitation de l'acide caprylique/du sulfate d'ammonium
<b>Taille:</b> 150ul , Concentration: 500 µg/ml by Bradford method using BSA as the standard;	<b>Identification du gène (NCBI):</b> 84280	<b>CloneNo.:</b> 4A1F8
<b>Hôte:</b> Mouse	<b>Nom complet:</b> BTB (POZ) domain containing 10	
<b>Isotype:</b> IgG2a	<b>MW calculé</b> 475 aa, 54 kDa	
<b>Immunogen Catalog Number:</b> AG10667	<b>MW observés:</b> 54 kDa	

## Applications

**Applications testées:**  
ELISA

**Spécificité de l'espèce:**  
Humain

## Informations générales

BTBD10 is a novel Akt-interacting protein with a BTB/POZ domain, which equally binds to any Akt. BTBD10 appears to behave as a suppressor of cell death including neuronal cell death related to ALS and an enhancer of cell growth via its positive regulation of Akt phosphorylation. The protein is expressed ubiquitously in vivo, distributing in cells as a unique cytoplasmic filamentous around the nucleus.

## Stockage

**Stockage:**  
Stocker à -20°C. Stable pendant un an après l'expédition.

**Tampon de stockage:**  
PBS avec azoture de sodium à 0,1 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20C

**\*\*\* Les 20ul contiennent 0,1% de BSA.**

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

