

À des fins de recherche uniquement

# Anticorps Monoclonal anti-FOXP3

Numéro de catalogue: 60149-1-Ig **1 Publications**



## Informations de base

Numéro de catalogue:

60149-1-Ig

Taille:

150ul , Concentration: 240 µg/ml by Bradford method using BSA as the standard;

Hôte:

Mouse

Isotype:

IgM

Numéro d'acquisition GenBank:

NM\_014009

Identification du gène (NCBI):

50943

Nom complet:

forkhead box P3

MW calculé

47 kDa

MW observés:

48 kDa

Méthode de purification:

Précipitation de l'acide caprylique/du sulfate d'ammonium

CloneNo.:

7B6A8

## Applications

Applications testées:

ELISA

Demandes citées:

IF

Spécificité de l'espèce:

Humain, rat, singe, souris

Espèces citées:

Humain

## Informations générales

FOXP3, also named as IPEX, JM2 and Scurfin, is a transcription factor. FOXP3 plays a critical role in the control of immune response. Defects in FOXP3 are the cause of immunodeficiency polyendocrinopathy, enteropathy, X-linked syndrome (IPEX); also known as X-linked autoimmunity-immunodeficiency syndrome. IPEX is characterized by neonatal onset insulin-dependent diabetes mellitus, infections, secretory diarrhea, thrombocytopenia, anemia and eczema. It is usually lethal in infancy. This antibody is a mouse monoclonal antibody raised against a peptide mapping within human FOXP3.

## Publications notables

Autrice	Pubmed ID	Journal	Application
Wakana Chikaishi	37220892	Thorac Cancer	IF

## Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20C

\*\*\* Les 20ul contiennent 0,1% de BSA.

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

