

À des fins de recherche uniquement

# Anticorps Monoclonal anti-LCK

Numéro de catalogue: 65028-1-Ig



## Informations de base

Numéro de catalogue:

65028-1-Ig

Taille:

150µL, 1333 µg/ml

Hôte:

Mouse

Isotype:

IgG2a

Numéro d'acquisition GenBank:

BC013200

Identification du gène (NCBI):

3932

Nom complet:

Lymphocyte-specific protein tyrosine kinase

MW calculé

539 aa, 56 kDa

MW observés:

50-60 kDa

Méthode de purification:

Purification par affinité contre l'antigène

Dilutions recommandées:

WB 1:500-1:2000

## Applications

Applications testées:

WB, ELISA

Spécificité de l'espèce:

Humain

Contrôles positifs:

WB : cellules Jurkat,

## Informations générales

LCK (Lymphocyte cell-specific protein-tyrosine kinase) is also named as LSK, p56-LCK and a nonreceptor protein-tyrosine kinase of the SRC oncogene family that is involved in transduction of T-cell receptor-mediated activation. It is implicated in the control of lymphocyte growth by virtue of its overexpression in some lymphoid malignancies and its transforming activity in heterologous systems (PMID:2784463). It has 3 isoforms produced by alternative splicing.

## Stockage

Stockage:

Stocker à une température comprise entre 2 °C et 8 °C.

Tampon de stockage:

PBS avec azoture de sodium à 0,1 % et glycérol à 50 % pH 7,3

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

