

À des fins de recherche uniquement

Anticorps Monoclonal anti-DR5

Numéro de catalogue: 65046-1-Ig



Informations de base

Numéro de catalogue:

65046-1-Ig

Taille:

150µl, 1000 µg/ml

Hôte:

Mouse

Isotype:

IgG1

Numéro d'acquisition GenBank:

BC001281

Identification du gène (NCBI):

8795

Nom complet:

tumor necrosis factor receptor
superfamily, member 10b

MW calculé

48 kDa

MW observés:

58 kDa

Méthode de purification:

Purification par protéine G

Applications

Applications testées:

ELISA

Spécificité de l'espèce:

Humain

Informations générales

TNFRSF10B, also known as CD262, DR5, TRAILR2, TRICK2 and KILLER, is a widely expressed single-pass type I membrane protein belonging to the tumour necrosis factor receptor superfamily (TNFRSF). It is a receptor for TNF-related apoptosis-inducing ligand (TRAIL), which is a member of the tumor necrosis factor (TNF) family of cytokines and induces apoptosis in a wide variety of cells (PMID: 9311998). TNFRSF10B contains two extracellular cysteine-rich repeats, typical for TNF receptor (TNFR) family members, and a cytoplasmic death domain (DD), through which TNFRSF10B is capable to transmit the apoptotic signal (PMID: 9311998; 20531300).

Stockage

Stockage:

Stocker à une température comprise entre 2 °C et 8 °C.

Tampon de stockage:

0.1M Tris-HCl, 0.15M NaCl, 0.04% merthiolate and 50% glycerol pH 7.4.

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

