

À des fins de recherche uniquement

Anticorps Monoclonal anti-CD24

Numéro de catalogue: 65121-1-Ig



Informations de base

Numéro de catalogue:

65121-1-Ig

Taille:

100ug , 500 µg/ml

Hôte:

Rat

Isotype:

IgG2b, kappa

Numéro d'acquisition GenBank:

BC075622

Identification du gène (NCBI):

12484

Nom complet:

CD24a antigen

Méthode de purification:

Purification par affinité

CloneNo.:

M1/69

Applications

Applications testées:

FC

Spécificité de l'espèce:

souris

Informations générales

Stockage

Stockage:

Store at 2-8°C. Stable for one year after shipment.

Tampon de stockage:

PBS avec azoture de sodium à 0,1 %, pH 7,3.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



1X10⁶ BALB/c mouse splenocytes were surface stained with 0.5 ug Anti-Mouse CD24 (65121-1-Ig, Clone: M1/69) and FITC anti-rat IgG2b Antibody at dilution 1:200 (green), or stained with 0.5 ug isotype control and FITC anti-rat IgG2b Antibody at dilution 1:200 (black). Cells were not fixed.