

À des fins de recherche uniquement

# Anticorps Monoclonal anti-CD1c

Numéro de catalogue: 65176-1-Ig



## Informations de base

Numéro de catalogue:

65176-1-Ig

Taille:

100ug, 1 mg/ml

Hôte:

Mouse

Isotype:

Mouse IgG1 kappa

Numéro d'acquisition GenBank:

BC126465

Identification du gène (NCBI):

911

Nom complet:

CD1c molecule

Méthode de purification:

Purifié par chromatographie d'affinité par protéine A

CloneNo.:

L161

## Applications

Applications testées:

FC

Spécificité de l'espèce:

Humain

## Informations générales

### Stockage

Stockage:

Stocker à une température comprise entre 2 °C et 8 °C.

Tampon de stockage:

PBS avec azoture de sodium à 0,1 %.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)


E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Données de validation sélectionnées



1X10<sup>6</sup> human peripheral blood lymphocytes were surface stained with APC Anti-Human CD19 (APC-65110, Clone: HIB19), and 1.5 ug Anti-Human CD1c (65176-1-Ig, Clone: L161) and FITC Anti-Mouse IgG1 Antibody. Cells were not fixed.



1X10<sup>6</sup> human peripheral blood lymphocytes were surface stained with APC Anti-Human CD19 (APC-65110, Clone: HIB19), and 1.5 ug Mouse IgG1 Isotype Control and FITC Anti-Mouse IgG1 Antibody. Cells were not fixed.