

À des fins de recherche uniquement

Anticorps Monoclonal anti-CD11a

Numéro de catalogue: 65177-1-Ig



Informations de base

Numéro de catalogue:

65177-1-Ig

Taille:

100ug, 0.5 mg/ml

Hôte:

Rat

Isotype:

IgG2a, kappa

Numéro d'acquisition GenBank:

BC145802

Identification du gène (NCBI):

16408

Nom complet:

integrin alpha L

Méthode de purification:

Purifié par chromatographie d'affinité par protéine A

CloneNo.:

M17/4

Applications

Applications testées:

FC

Spécificité de l'espèce:

souris

Informations générales

CD11a, also known as integrin α L or LFA-1 α , is a 170-180 kDa transmembrane glycoprotein that forms a heterodimer with CD18 (PMID: 7027264; 1672643). CD11a/CD18 (LFA-1) is expressed by all leukocytes and mediates cell adhesion through interactions with its ligands, intercellular adhesion molecule 1 (ICAM-1), ICAM-2, and ICAM-3 (PMID: 7479767). CD11a/CD18 also functions in lymphocyte costimulatory signaling (PMID: 1972160).

Stockage

Stockage:

Store at 2-8°C. Stable for one year after shipment.

Tampon de stockage:

PBS with 0.09% sodium azide.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



1X10⁶ mouse thymocytes were surface stained with 0.5 ug Anti-Mouse CD11a (65177-1-Ig, Clone:M17/4) (red), or 0.5 ug Isotype Control. Cells were not fixed.