

À des fins de recherche uniquement

Anticorps Monoclonal anti-CD41

Numéro de catalogue: 65189-1-Ig



Informations de base

Numéro de catalogue:

65189-1-Ig

Taille:

100ug, 0.5 mg/ml

Hôte:

Rat

Isotype:

IgG1, kappa

Numéro d'acquisition GenBank:

BC120493

Identification du gène (NCBI):

16399

Nom complet:

integrin alpha 2b

Méthode de purification:

Purifié par chromatographie d'affinité par protéine G

CloneNo.:

MWReg30

Applications

Applications testées:

FC

Spécificité de l'espèce:

souris

Informations générales

Stockage

Stockage:

Store at 2-8°C. Stable for one year after shipment.

Tampon de stockage:

PBS with 0.09% sodium azide.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



1X10⁶ mouse splenocytes were surface stained with 0.5 ug Anti-Mouse CD41 (65189-1-Ig, Clone:MWReg30) (red), or 0.5 ug Isotype Control. Cells were not fixed.