

À des fins de recherche uniquement

Anticorps Monoclonal anti-CD4

Numéro de catalogue: 65199-1-Ig



Informations de base

Numéro de catalogue:

65199-1-Ig

Taille:

500ug, 0.5 mg/ml

Hôte:

Rat

Isotype:

IgG2b, kappa

Numéro d'acquisition GenBank:

BC039137

Identification du gène (NCBI):

12504

Nom complet:

CD4 antigen

Méthode de purification:

Purification par protéine A

CloneNo.:

RM4-4

Applications

Applications testées:

FC

Spécificité de l'espèce:

souris

Informations générales

CD4 is a 55-kDa transmembrane glycoprotein expressed on T helper cells, majority of thymocytes, monocytes, macrophages, and dendritic cells (PMID: 9304802; 12213222). CD4 is an accessory protein for MHC class-II antigen/T-cell receptor interaction. It plays an important role in T helper cell development and activation (PMID: 9539765; 3112582). CD4 serves as a receptor for the human immunodeficiency virus (HIV) (PMID: 9304802).

Stockage

Stockage:

Store at 2-8°C. Stable for one year after shipment.

Tampon de stockage:

PBS with 0.09% sodium azide.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



1X10⁶ mouse splenocytes were surface stained with 0.25 ug Anti-Mouse CD4 (65199-1-Ig, Clone:RM4-4) and FITC anti-rat IgG2b Antibody at dilution 1:100 (red), or 0.25 ug Control Antibody. Cells were not fixed.