

À des fins de recherche uniquement

Anticorps Monoclonal anti-IgG1 Isotype Control



Numéro de catalogue: 65212-1-Ig

Informations de base

Numéro de catalogue:

65212-1-Ig

Taille:

500ug, 0.5 mg/mL

Hôte:

Rat

Isotype:

IgG1, κ

Numéro d'acquisition GenBank:

Identification du gène (NCBI):

Nom complet:

Méthode de purification:

Purification par affinité

CloneNo.:

HRPN

Applications

Applications testées:

FC

Spécificité de l'espèce:

Informations générales

The HRPN immunoglobulin is useful as an isotype-matched control. The HRPN immunoglobulin has an unknown binding specificity and is used as an isotype control for rat IgG1 antibodies. This antibody has been quality-tested for flow cytometry as negative control.

Stockage

Stockage:

Stocker à une température comprise entre 2 °C et 8 °C.

Tampon de stockage:

Tampon à base de phosphate avec azoture de sodium à 0,09 %, pH 7,2.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



1X10⁶ mouse splenocytes were surface stained with 0.5 μ g Rat IgG1 Isotype Control (65212-1-Ig, Clone: HRPN) and FITC anti-rat IgG1 Antibody at dilution 1:200 (green), or not stained (blank) (black). Cells were not fixed.