

À des fins de recherche uniquement

# Anticorps Monoclonal anti-CD48

Numéro de catalogue: 65227-1-Ig



## Informations de base

Numéro de catalogue:

65227-1-Ig

Taille:

500ug, 0.5 mg/mL

Hôte:

Armenian Hamster

Isotype:

IgG

Numéro d'acquisition GenBank:

BC060977

Identification du gène (NCBI):

12506

Nom complet:

CD48 antigen

Méthode de purification:

Purification par affinité

CloneNo.:

HM48-1

## Applications

Applications testées:

FC

Spécificité de l'espèce:

souris

## Informations générales

### Stockage

Stockage:

Store at 2-8°C. Stable for one year after shipment.

Tampon de stockage:

PBS with 0.09% sodium azide.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Données de validation sélectionnées



1X10<sup>6</sup> mouse splenocytes were surface stained with 0.5 ug Anti-Mouse CD48 (65227-1-Ig, Clone:HM48-1) and FITC anti-Armenian Hamster IgG Antibody at dilution 1:100 (red), or 0.5 ug Control Antibody. Cells were not fixed.