

À des fins de recherche uniquement

Anticorps Monoclonal anti-CD55

Numéro de catalogue: 65271-1-Ig



Informations de base

Numéro de catalogue:

65271-1-Ig

Taille:

100ug, 0.5 mg/ml

Hôte:

Mouse

Isotype:

IgG1

Numéro d'acquisition GenBank:

BC001288

Identification du gène (NCBI):

1604

Nom complet:

CD55 molecule, decay accelerating factor for complement (Cromer blood group)

MW calculé

41 kDa

Méthode de purification:

Purification par affinité

CloneNo.:

JS11

Applications

Applications testées:

FC

Spécificité de l'espèce:

Humain

Informations générales

Stockage

Stockage:

Store at 2-8°C. Stable for one year after shipment.

Tampon de stockage:

PBS with 0.09% sodium azide.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)

E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



1X10⁶ human PBMCs were surface stained with 0.2 ug Anti-Human CD55 (65271-1-Ig, Clone:JS11K) and CoraLite®594-Conjugated AffiniPure Goat Anti-Mouse IgG(H+L) at dilution 1:1000. Cells were not fixed. Lymphocytes were gated.