

À des fins de recherche uniquement

Anticorps Monoclonal anti-RXRG

Numéro de catalogue: 66306-1-Ig **1 Publications**



Informations de base

Numéro de catalogue: 66306-1-Ig	Numéro d'acquisition GenBank: BC012063	Méthode de purification: Purification par protéine A
Taille: 150ul, Concentration: 1300 µg/ml by Nanodrop and 1000 µg/ml by Bradford method using BSA as the standard;	Identification du gène (NCBI): 6258	CloneNo.: 1E7E5
Hôte: Mouse	Nom complet: retinoid X receptor, gamma	Dilutions recommandées: WB 1:1000-1:4000
Isotype: IgG2b	MW calculé: 51 kDa	
Immunogen Catalog Number: AG1611	MW observés: 51 kDa	

Applications

Applications testées: WB, ELISA	Contrôles positifs: WB : cellules A549,
Spécificité de l'espèce: Humain	
Espèces citées: Humain	

Informations générales

RXRG, also named as, Retinoic acid receptor RXR-gamma, is a 463 amino acid protein, which localize in the nucleus and may form homodimer. Retinoic acid receptors bind as heterodimers to their target response elements in response to their ligands, all-trans or 9-cis retinoic acid, and regulate gene expression in various biological processes. The RAR/RXR heterodimers bind to the retinoic acid response elements (RARE) composed of tandem 5'-AGGTCA-3' sites known as DR1-DR5.

Publications notables

Autrice	Pubmed ID	Journal	Application
Qian Hao	33627621	Cell Death Dis	

Stockage


Stockage:
Stocker à -20°C. Stable pendant un an après l'expédition.
Tampon de stockage:
PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3
L'aliquotage n'est pas nécessaire pour le stockage à -20C

*** Les 20ul contiennent 0,1% de BSA.

For technical support and original validation data for this product please contact:
T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)
E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Données de validation sélectionnées



A549 cells were subjected to SDS PAGE followed by western blot with 66306-1-Ig (RXRG Antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.