

À des fins de recherche uniquement

# Anticorps Monoclonal anti-SIX2

Numéro de catalogue: 66347-1-Ig **5 Publications**



## Informations de base

Numéro de catalogue: 66347-1-Ig	Numéro d'acquisition GenBank: BC024033	Méthode de purification: Purification par protéine A
Taille: 150ul, Concentration: 1600 µg/ml by Nanodrop and 1000 µg/ml by Bradford method using BSA as the standard;	Identification du gène (NCBI): 10736	CloneNo.: 1F9G11
Hôte: Mouse	Nom complet: SIX homeobox 2	Dilutions recommandées: WB 1:1000-1:6000 IHC 1:50-1:500
Isotype: IgG2b	MW calculé: 291 aa, 32 kDa	
Immunogen Catalog Number: AG2124	MW observés: 37 kDa	

## Applications

### Applications testées:

IHC, WB, ELISA

### Demandes citées:

IF, IHC, WB

### Spécificité de l'espèce:

canin, Humain, porc, rat, souris

### Espèces citées:

souris

**Remarque-IHC: il est suggéré de démasquer l'antigène avec un tampon de TE buffer pH 9,0; (\*) À défaut, 'le démasquage de l'antigène peut être 'effectué avec un tampon citrate pH 6,0.**

### Contrôles positifs:

WB : tissu de muscle squelettique de porc, cellules HEK-293, cellules LO2, cellules MDCK, cellules NCCIT, tissu de muscle squelettique de rat, tissu de muscle squelettique de souris

IHC : tissu embryonnaire de souris,

## Informations générales

The SIX proteins (sine oculis) are a family of homeodomain transcription factors that share a conserved DNA binding domain. Six members (Six1-Six6) of the Six gene family have been identified in mice and humans. SIX2, containing one homeobox DNA-binding domain, is highly expressed in fetal tissues but expression is limited in adult tissues. SIX2 may be involved in limb tendon and ligament development[PMID:21420949]. It has been previously shown that SIX2 is expressed in developing mesenchymal tissue including head and urogenital system at the time of overt midfacial and renal differentiation[PMID: 22282599].

## Publications notables

Autrice	Pubmed ID	Journal	Application
Jun Li	36130284	Nucleic Acids Res	IHC
Wen Shi	36057966	Dev Dyn	IHC
Bliss Magella	29183737	Dev Biol	IHC,IF

## Stockage

### Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

### Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20C

\*\*\* Les 20ul contiennent 0,1% de BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com  
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


## Données de validation sélectionnées



Various cells were subjected to SDS PAGE followed by western blot with 66347-1-Ig (SIX2 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.



Various tissues were subjected to SDS PAGE followed by western blot with 66347-1-Ig (SIX2 antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.



Immunohistochemical analysis of paraffin-embedded mouse embryo tissue slide using 66347-1-Ig (SIX2 antibody) at dilution of 1:200 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).