

À des fins de recherche uniquement

# Anticorps Monoclonal anti-THRA

Numéro de catalogue: 66703-1-Ig **Phare**



## Informations de base

Numéro de catalogue:

66703-1-Ig

Taille:

150ul

Hôte:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG27189

Numéro d'acquisition GenBank:

BC002728

Identification du gène (NCBI):

7067

Nom complet:

thyroid hormone receptor, alpha  
(erythroblastic leukemia viral (v-erb-  
a) oncogene homolog, avian)

MW calculé

54 aa, 7 kDa

Méthode de purification:

Purification par protéine G

CloneNo.:

1A11C8

## Applications

Applications testées:

ELISA

Spécificité de l'espèce:

Humain, rat, souris

## Informations générales

### Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azoture de sodium à 0,1 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20C

\*\*\* Les 20ul contiennent 0,1% de BSA.

For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free  
in USA), or 1(312) 455-8498 (outside USA)

E: [proteintech@ptglab.com](mailto:proteintech@ptglab.com)  
W: [ptglab.com](http://ptglab.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Données de validation sélectionnées



WB result of THRA antibody (66703-1-Ig; 1:4000;  
incubated at room temperature for 1.5 hours) with  
sh-Control and sh-THRA transfected A431 cells.