

À des fins de recherche uniquement

Anticorps Monoclonal anti-STRAP

Numéro de catalogue: 66712-1-Ig

Phare

1 Publications



Informations de base

Numéro de catalogue: 66712-1-Ig	Numéro d'acquisition GenBank: BC000162	Méthode de purification: Purification par protéine G
Taille: 150ul, Concentration: 1300 µg/ml by Nanodrop and 1000 µg/ml by Bradford method using BSA as the standard;	Identification du gène (NCBI): 11171	CloneNo.: 5E11F3
Hôte: Mouse	Nom complet: serine/threonine kinase receptor associated protein	Dilutions recommandées: WB 1:5000-1:20000 IHC 1:250-1:1000 IF 1:50-1:500
Isotype: IgG1	MW calculé: 38 kDa	
Immunogen Catalog Number: AG13085	MW observés: 38 kDa	

Applications

Applications testées:

IF, IHC, WB, ELISA

Demandes citées:

WB

Spécificité de l'espèce:

Humain, rat, souris

Espèces citées:

Humain

Remarque-IHC: il est suggéré de démasquer l'antigène avec un tampon de TE buffer pH 9,0; (*) A défaut, 'le démasquage de l'antigène peut être effectué avec un tampon citrate pH 6,0.

Contrôles positifs:

WB : cellules HeLa, cellules HEK-293, cellules HepG2, cellules HSC-T6, cellules Jurkat, cellules NIH/3T3

IHC : tissu ovarien de souris,

IF : cellules HEK-293,

Informations générales

STRAP (Serine-threonine kinase receptor-associated protein), also known as MAWD, is a ubiquitous WD40 domain protein, common function of which is to provide a suitable scaffold for coordinating multiprotein complex assemblies and thus, regulate a variety of cellular processes like signal transduction, transcriptional regulation, programmed cell death and so on. STRAP plays a role in the cellular distribution of the SMN complex, which is essential for spliceosomal snRNP assembly in the cytoplasm. STRAP inhibits transforming growth factor-beta (TGF-beta) signaling and enhances tumorigenicity via TGF-beta-dependent and -independent mechanisms. STRAP is up-regulated in several cancers and functions as an oncogene. STRAP imparts oncogenic characteristics to cells by promoting ERK and pRb phosphorylation. Moreover, STRAP regulates c-Jun stability by decreasing the ubiquitylation and proteosomal degradation of c-Jun.

Publications notables

Autrice	Pubmed ID	Journal	Application
Tian-Tian Xu	34504787	Front Oncol	WB

Stockage

Stockage:

Stocker à -20°C. Stable pendant un an après l'expédition.

Tampon de stockage:

PBS avec azoture de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20C

*** Les 20ul contiennent 0,1% de BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Données de validation sélectionnées




Various lysates were subjected to SDS PAGE followed by western blot with 66712-1-Ig (STRAP antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.




WB result of STRAP antibody (66712-1-Ig; 1:5000; incubated at room temperature for 1.5 hours) with sh-Control and sh-STRAP transfected HEK-293 cells.



Immunohistochemical analysis of paraffin-embedded mouse ovary tissue slide using 66712-1-Ig (STRAP antibody) at dilution of 1:500 (under 10x lens. Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunohistochemical analysis of paraffin-embedded mouse ovary tissue slide using 66712-1-Ig (STRAP antibody) at dilution of 1:500 (under 40x lens. Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunofluorescent analysis of (4% PFA) fixed HEK-293 cells using STRAP antibody (66712-1-Ig, Clone: 5E11F3) at dilution of 1:200 and Coralite@488-Conjugated AffiniPure Goat Anti-Mouse IgG(H+L).