

À des fins de recherche uniquement

Anticorps Monoclonal anti-PRTN3

Numéro de catalogue: 67030-1-Ig



Informations de base

Numéro de catalogue: 67030-1-Ig	Numéro d'acquisition GenBank: BC096183	Méthode de purification: Purification par protéine A
Taille: 150ul, Concentration: 2100 µg/ml by Nanodrop and 1000 µg/ml by Bradford method using BSA as the standard;	Identification du gène (NCBI): 5657	CloneNo.: 2D1B8
Hôte: Mouse	Nom complet: proteïnase 3	Dilutions recommandées: WB 1:1000-1:4000 IHC 1:1000-1:4000 IF 1:200-1:800
Isotype: IgG2b	MW calculé: 256 aa, 28 kDa	
Immunogen Catalog Number: AG25915	MW observés: 28 kDa	

Applications

Applications testées:

IF, IHC, WB, ELISA

Spécificité de l'espèce:

Humain

Remarque-IHC: il est suggéré de démasquer l'antigène avec un tampon de TE buffer pH 9.0; (*) À défaut, 'le démasquage de l'antigène peut être effectué avec un tampon citrate pH 6,0.

Contrôles positifs:

WB : tissu placentaire humain,

IHC : tissu splénique humain, tissu d'amygdales humaine

IF : tissu d'amygdales humaine,

Informations générales

Stockage

Stockage:

Stocker à -20 °C.

Tampon de stockage:

PBS avec azote de sodium à 0,02 % et glycérol à 50 % pH 7,3

L'aliquotage n'est pas nécessaire pour le stockage à -20C

*** Les 20ul contiennent 0,1% de BSA.


For technical support and original validation data for this product please contact:

T: 1 (888) 4PTGLAB (1-888-478-4522) (toll free
in USA), or 1(312) 455-8498 (outside USA)


E: proteintech@ptglab.com
W: ptglab.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Données de validation sélectionnées




human placenta tissue were subjected to SDS PAGE followed by western blot with 67030-1-Ig (PRTN3 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.



Immunohistochemical analysis of paraffin-embedded human spleen tissue slide using 67030-1-Ig (PRTN3 antibody) at dilution of 1:2000 (under 10x lens. Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunohistochemical analysis of paraffin-embedded human spleen tissue slide using 67030-1-Ig (PRTN3 antibody) at dilution of 1:2000 (under 40x lens. Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunofluorescent analysis of (4% PFA) fixed human tonsillitis tissue using PRTN3 antibody (67030-1-Ig, Clone: 2D1B8) at dilution of 1:400 and CoraLite®488-Conjugated AffiniPure Goat Anti-Mouse IgG(H+L).